



植絨纖維拭子

Product Reference	Product Description	Length of Swab Shaft	Diameter of Swab Shaft	Immediate packaging	Sterility
MW002NF	Dryswab™, Nylon Flock Fibre, Mini-Tip & Flexible Shaft, Breakpoint 78mm	151mm with 78mm breakpoint	2.5mm (Tapers to 1.0mm at tip)	Individual in Peel Pouch	Sterilised by Gamma Irradiation

Product Description

Product code MW002NF 植絨纖維 拭子

在 78 毫米處有一個斷點，2.5mm（逐漸變細到1.0mm 尖端）單支包裝滅菌伽瑪輻照滅菌與許多拭子採集瓶兼容，




Figure 1 Peel pouch example 1





Figure 2 Peel pouch example 2



MW002NF

常規QC測試規範

每一批都經過測試，以確保流感嗜血桿菌 NCTC 8143在指定的持續時間恢復生長沒有被抑制。Dryswab® 產品無法通過 QC 性能測試時間，取決於被測試的產品和微生物，經過指定的孵育時間後與 0 小時的對照管相比結果下降了，超過 2-log。

監管

所有 Dryswab™ 均帶有 CE 標誌並符合：歐洲醫療器械指令(93/42/EEC)，1993，歐洲共同體官方期刊，L169

質量保證

醫療電線和設備的設計已通過 ISO 9001:2015、ISO 13485:2016 認證無菌醫療器械的設計和製造，以及無菌和非無菌體外診斷設備。