























Special Application

			Bud	Shaft	Flexible	Break Point	Length	Pack Size	
Cytology		Brush swab for specimen collection in cytology. The growing importance of routine screening for cervical cancer has required the development of improved methods for cell collection. Cytotak™ is a swab-length cytology collection brush which will give excellent results in the collection of endocervical cells. The tapered spiral configuration assists access to narrow opening and gentle rotation creates a brushing action against skin surface for enhanced recovery of cellular material. The product is CE-marked, and is available in three convenient formats.							
	MW146	Cytotak™ cytology sampling brush, bulk, non-sterile	Brush	Plastic	No	No	Standard	100	
	MW147	Cytotak™ cytology sampling brush, in Peel Pouch	Brush	Plastic	No	No	Standard	125	
	MW148	Cytotak™ cytology sampling brush, in labelled tube	Brush	Plastic	No	No	Standard	100	
Podiatry	Taper-tipped swabs for under-nail sampling.								
	MW109	Dryswab, Taper tipped, in Peel Pouch	Cotton	Wood	No	No	Standard	1250	

Veterinary		The range of Dry Swabs includes several longer swabs particularly suited to larger animals. In some of these the long shaft is enclosed in a silicone sheath into which the swab head can be with withdrawn as the swab is retrieved from the sampling site.							
	MW123	Equine Flu Swab, 40cm, in Pouch	Rayon	Plastic	Yes	No	40cm	100	
	MW123/32	Equine Flu Swab, 32cm, in Pouch	Rayon	Plastic	Yes	No	32cm	100	
	MW124	Nasopharyngeal Swab with bud for large animal	Rayon	Plastic	Yes	No	27cm	100	
	MW126	Nasopharyngeal Swab with brush for large animal	Brush	Twisted wire	Yes	No	27cm	100	
	MW128	Laryngeal Swab	Rayon	Straight wire	Yes	No	27cm	100	
	MW100	Dryswab™, Fine Tip	Rayon	Plastic	Yes	No	Standard	125	
	MW113	Dryswab™, Fine Plastic Stick	Rayon	Plastic	Yes	No	Standard	125	
	MW142	Dryswab™ ENT, Straight Wire	Rayon	Straight wire	No	No	Standard	100	

Bulk Packed

			Bud	Shaft	Flexible	Break Point	Length	Pack Size	
PurFlock	MW810B	PurFlock® Ultra, Standard Tip, Tapered Shaft	PurFlock®	Plastic	No	80mm	Standard	1000	
	MW811B	PurFlock® Ultra, Standard Tip, Flexible Shaft	PurFlock®	Plastic	Yes	80mm	Standard	1000	
	MW812B	PurFlock® Ultra, Mini Tip, Flexible Shaft	PurFlock®	Plastic	Yes	No	Standard	1000	
	MW813B	PurFlock® Ultra, Ultrafine Tip, Flexible Shaft	PurFlock®	Plastic	Yes	100mm	Standard	1000	
	MW814B	PurFlock® Ultra, Micro Ultrafine Tip, Flexible Shaft	PurFlock®	Plastic	Yes	100mm	Standard	1000	
HydraFlock	MW816B	HydraFlock®, Standard Tip, Tapered Shaft	HydraFlock®	Plastic	No	80mm	Standard	1000	
	MW817B	HydraFlock®, Standard Tip, Flexible Shaft	HydraFlock®	Plastic	Yes	80mm	Standard	1000	
	MW818B	HydraFlock®, Mini Tip, Flexible Shaft	HydraFlock®	Plastic	Yes	No	Standard	1000	
	MW819B	HydraFlock®, Ultrafine Tip, Flexible Shaft	HydraFlock®	Plastic	Yes	100mm	Standard	1000	
	MW820B	HydraFlock®, Micro Ultrafine Tip, Flexible Shaft	HydraFlock®	Plastic	Yes	100mm	Standard	1000	
Speciality	MW121	Dryswab™, Fine Plastic Stick	Rayon	Plastic	No	No	Standard	1000	
	MW131	Dryswab™ ENT, Straight wire	Rayon	Straight wire	No	No	Standard	1000	
	MWSTICKS	Wooden applicator sticks, Non-Sterile	No bud	Wood	No	No	Standard	5000	