



FastGene™



FastGene™

NANOVIEW PHOTOMETER

Microvolume & Cuvette

DISCOVER NEW HORIZONS

1

TWO SOURCES MICROVOLUME AND CUVETTE READER

Drop or cuvette - Why not both in one device?
The FastGene® NanoView comes with a microvolume drop reader and a cuvette reader. For easy analysis of small and large sample volumes.

2

RELIABLE ANALYSIS THE SPECIALIST

Three distinct measurement wavelengths (260 nm, 280 nm, 600 nm) enable precise measurements of DNA/RNA, proteins and bacterial cultures.

3

INTUITIVE USAGE EASY TO USE SOFTWARE

The software of the FastGene® NanoView was made for easy handling. The touchscreen enables quick and intuitive navigation through different measurement modes.



ABOUT US

NIPPON Genetics is a Japanese-German Life-Tech company, which focuses on the development of cutting edge products. We have over 30 years experience in the life-tech sector with offices in Tokyo and Kyoto, Japan, as well as in Dueren, Germany.

Need an online product demonstration?

We understand that, as a scientist, you would like to test the FastGene® NanoView before buying it. That is why we offer a product demonstration online or in your lab.

Just arrange an appointment with us!




Visit our website for other products, such as our MIDORI^{Green} DNA/RNA dyes, Nucleic Acid Purification kits, PCR enzymes, Instruments and many more...



www.nippongenetics.eu

**Drop or cuvette -
 Why not both in one
 device?**

CUVETTE READER



MICROVOLUME READER



Cat. No.: FG-NP02

**NEED AN ONLINE PRODUCT
 DEMONSTRATION?**

We understand that, as a scientist,
 you would like to test the
 FastGene® NanoView before buying it.

That is why we offer a product
 demonstration online or in your lab.


Just arrange an appointment
 with us!



THE COMPACT SPECIALIST PHOTOMETER

- ✔ **Best device performance** - Microvolume drop and cuvette reader
- ✔ **Specialist** - For DNA, RNA, Proteins and OD600
- ✔ **Convenient usage** - 10 preset measurement modes
- ✔ **Smart** - Compact, silent and easy to use
- ✔ **Cost-saving** - Affordable stand-alone device

PRECISE



RELIABLE SAMPLE ANALYZER

The specialist device comes with 10 preset measurement modes enabling decisive DNA, RNA, protein or OD600 analysis.


CONVENIENT



SIMPLE TO USE SOFTWARE

The software UI of the FastGene® NanoView was designed for ease of use with large and self-explanatory icons.

SMART



STAND-ALONE DEVICE

The microvolume drop reader and cuvette reader are both integrated under one lid making the FastGene® NanoView ultra-small, silent and an absolute lightweight.

OPTIMIZED INTERFACE - EASY MEASUREMENT



The software of the FastGene® NanoView was made for easy handling. The main menu tab enables easy navigation through distinct measurement modes via touchscreen with intuitive icons.